

# VESSELS AND PROCESSING SYSTEMS



- PRESSURE VESSELS FROM 0.2 TO 1000 LITERS
- COMPLETE CUSTOM TURNKEY PROCESSING AND LIQUID HANDLING SYSTEMS DESIGNED TO SUIT CLIENT PROCESSING APPLICATION REQUIREMENTS
- INTEGRATED WITH VALVING, MIXERS, SENSORS, HEATERS, CONDENSERS, PUMPS, INSTRUMENTATION AND CONTROL SYSTEMS, ETC.
- SPECIAL DESIGNS FOR REACTORS, FERMENTERS, MIXERS, HOMOGENIZERS, DISPENSERS, HYDROGENATORS, ETC.
- HIGH PURITY CHEMICAL, PHARMACEUTICAL, FOOD USE FABRICATION AND FINISHES
- STAINLESS STEEL, HASTELLOY AND OTHER MATERIALS, ASME STANDARD, TSSA, CE/PED, UN/DOT, ETC.

## NUTSCHE FILTER-DRYER SYSTEMS

- PURIFICATION AND COLLECTION OF HIGH PURITY SOLIDS, PHARMACEUTICAL INTERMEDIATES, ETC.
- SINGLE VESSEL FILTERING, WASHING, MIXING, RESLURRYING AND DRYING; 0.2 TO 500L SIZES
- A LOGICAL NEXT STEP UP FROM BUCHNER FUNNEL LABWORK, WITH PRESSURIZATION, VACUUM, AGITATION, HEATING
- STAINLESS STEEL, HASTELLOY, OTHER MATERIALS, ASME CERTIFIED, SANITARY SURFACES
- WIDE RANGE OF DESIGNS, MIXERS, AGITATORS, MOUNTINGS, FILTER TYPES, INSTRUMENTATION, HARVESTING, ETC.
- NEW STANDARD *BENCHTOP NUTSCHE FILTER SYSTEMS*



**Solution Driven.**