

Pope Laboratory Scale Wiped-Film Stills



Pope 2" Wiped-Film Molecular Still with typical accessories:
External Condenser with Receiver, Vacuum Gauge, RPM Indicator, Metal Band Heater (insulation removed),
Temperature Controller, and Vacuum Pump System (high-vacuum type with oil diffusion booster)

Distillation and Evaporation of Heat Sensitive Materials

2", 4" AND 6" DIAMETER MODELS

Diam. inches	Surface Area ft ² (m ²)	Typical Flowrates Lit/hr
2	0.35 (0.033)	0.1-3
4	1.16 (0.108)	1-20
6	2.35 (0.218)	2-40
6	3.53 (0.327)	3-60

Pope Scientific, Inc.
351 N. Dekora Woods Blvd.
P.O. Box 80018
Saukville, WI 53080 USA
Phone: 262-268-9300 • Fax: 262-268-9400
www.popeinc.com • Dewars – dewars@popeinc.com
Corporate Capabilities and others – info@popeinc.com

